



SPECTRUM - MS 400 Electrosurgical Unit



ELECTRO-SURGICAL UNIT

Main features

- Widely applicable for ordinary surgery, cardiology, urology, gynecology, underwater cut
- Automatic performance check during startup
- System that provides surgeons with improved performance at lower power settings, minimizing the risk of tissue damage
- Max 400W power output
- 6 working mode: pure/blend/cautery cut, spray/contact coagulation, bipolar coagulation
- Power output automatically cut off in any case of interruption
- Operated by hand or foot pedal

OPTIONS

- BIPOLAR FORCEP
- TROLLY

FRONT VIEW



REAR VIEW

