

Fetal Beat S-12

(7 inch Multi-color LCD)

Main Features:

- Vertical Type Design
- 7 inch vivid Display
- Single or Twin FHR (Ultrasound) Transducers
- Data Graph and Trend table review
- Thermal Printer with 110mm wide paper
- Central monitoring system (optional)
- Programmable Alarms



Specification:

FHR:

Frequency: 1MHZ
Range: 30- 240 BPM
Resolution: 1BPM
Accuracy: ± 2 BPM

Printer:

Thermal Printer with high print density (1000dpi)
Recording Speed 10, 20, 30 mm/min
Recording Paper: 110mm roll

TOCO:

Range: 0 - 100 %
Resolution: 1
Non-linear Error: $\leq \pm 8\%$
Zero Control: Manual

Physical Specification:

Height: 28.7cm
Width: 18.6 cm
Length: 22.0 cm

Operation Environment:

Temperature: 10 °C - 40 °C
Relative Humidity: $\leq 70\%$

Power Supply:

AC 100 - 240V 60Hz
Built-in rechargeable battery