

OXICARE VS 2000 SN

3 PARAMETERS PR + SpO2 + NIBP

TABLE-TOP OXIMETER

can continually measure SpO2, PR at a real time. It is with compact appearance, comprehensive functions, simple and convenient operation, and widely used in hospitals, community medical treatment and family.

VS 2000 S

2 PARAMETERS PR + SpO2



Features:

- High resolution 4.3" color LCD screen
- Compact and portable, allowing for uninterrupted monitoring
- AC power supply, with built-in lithium battery
- Adjustable audible and visual alarms
- Storage and review of patient's trend data up to 36 hours
- Applicable to adult, pediatric and neonatal patient
- Convenient and quick in inquiring measurement data
- USB data upload
- Optional: Printing

Performance:

- NIBP Method: Oscillometry
- Operation modes: Manual/Automatic/STAT
- Auto Measure time : 5~250 minutes adjustable, interval 5 minutes
- Measurement range: 10~270mmHg
- Overpressure protection: yes
- Resolution: 1mmHg
- Alarm: systolic, Diastolic, Mean
- SPO2 Measurement range: 0~100%
- Resolution: 1%
- Accuracy: 70~100% $\pm 2\%$, 0~69% unspecified
- Alarm range: 0~100%
- Pulse Rate Range: 0~254bpm
- Resolution: 1bpm
- Accuracy: ± 1 bpm